



BIOFIX

ADHESIVE FOR BRACKETS





- **Adhesive for brackets**
- **Light Cured**
- **Ideal consistence for the brackets fixation, does not allow the bracket to slide after positioned**
- **Presents fluorine guaranteeing a prophylactic action against demineralization in the bracket area**





INDICATION

- For fixation of metallic, ceramic and/or plastic orthodontic bracket in the dental surface



PRESENTATION

- 2 x 2,5g syringes of BIOFIX
- 1 syringe of ATTAQUE GEL
- 5 application tips





BENEFITS

- **Light Cured**
- **Contain fluorine**
- **Great consistency**
- **Excellent adhesion to the dental enamel**
- **Easy removal and application**

