




INTERIM

Temporary Restorative Material





- **Temporary restorative material, based on Zinc Oxide and Eugenol**
- **NEW formulation, easier to handle**
- **Setting Time: from 5 to 8 minutes.**
- **The powder is present in ivory color**
- **Resistant to compression**
- **High marginal sealing**
- **Antiseptics and sedatives properties due to Eugenol**
- **Insoluble in the mouth**





INDICATION

- **Indicated for cavities lining and long time restorations**





PRESENTATION

- 1x38g flask of INTERIM powder
- 1x15ml flask of INTERIM liquid



Powder



Liquid



BENEFITS

- **Easy to handle, easy to insert in the cavity**
- **Low setting time, more agility**
- **Low marginal infiltrations indices**
- **Low postoperative pain indices**
- **Biocompatibility with pulp protector**

